


## **Quick Start Manual**

3

2.







Re-mount Cooling Unit (with Mouthpiece).

Switch on Vaporizer. Ready to use when actual temperature has reached set temperature.

Inhale from the Mouthpiece.

Refer to the Instructions for Use for important safety and user advice. Do not operate the device without first reading and understanding the Instructions for Use!

STORZ & BICKEL GmbH · In Grubenäcker 5-9 · 78532 Tuttlingen/Germany



## Dismantling, Cleaning, Replacing





Cooling Unit





Pull the Cap Lock towards the Mouthpiece.





Pull the Cooling Unit Cap off the Vaporizer.

7.







3.









Pull the Mouthpiece out of the Cooling Unit Cap.

Unscrew the Cooling Unit Base from the Vaporizer.

Lever out the large Base Seal Ring with the Filling Chamber Tool.

Press out the Screen.

Refer to the Instructions for Use for important safety and user advice. Do not operate the device without first reading and understanding the Instructions for Use!

STORZ & BICKEL GmbH · In Grubenäcker 5-9 · 78532 Tuttlingen/Germany